

CFPM400 & CFPM401

4 Movement Tables

- **Industry's Widest Adjustable Height Range 26"- 44"**
- 56" Metal Free Imaging Space
- Cost Effective 4 Movement Table
- 500 lbs. Patient Weight Capacity
- **Warranty: 4 years parts, 1 year labor**



RECTANGULAR TOP OPTION



See OAKWORKS®' CFPM400 in action!

Use your smartphone to view the CFPM400 video or watch online at www.oakworksmed.com.

CFPM400 & CFPM401 PRODUCT SPECS

PATIENT WEIGHT CAPACITY	500 lbs. (227kg)
CARBON FIBER TOP	Integrated Headrest 22" x 84" (56cm x 214cm) Rectangular 24" x 80" (61cm x 203cm)
CFPM400 - MOVEMENTS	
MOTION 1	height range 26" - 44" (66cm - 112cm)
MOTION 2	Trendelenburg tilt +15°/-12°
MOTION 3	10" (25cm) longitudinal travel
MOTION 4	lateral tilt ±15°
CFPM401 - MOVEMENTS	
MOTION 1	height range 26" - 44" (66cm - 112cm)
MOTION 2	Trendelenburg tilt +15°/-11°
MOTION 3	10" (25cm) longitudinal travel
MOTION 4	±4" (10cm) lateral travel
CONTROLS	1 foot (4 motion), 1 hand (4 motion)
TABLE WEIGHT	475 lbs. (216 kg.)
WARRANTY	4 years parts, 1 year labor
CERTIFICATIONS	FDA listed, cETL approved, CE marked



Certifications: UL 60601-1; CSA C22.2 # 601.1, IEC 60601-1; IEC 60601-2-32 IEC 60601-1-2

